

Data Sheet

Water for Injection (WFI) / Water for Cell Culture (WFI) WFI Quality Water

Product Name:

1A Clean Water for Injection (WFI)

Size: 200 Liter 3D Bag

500 Liter 3D Bag

Product Name:

1A Clean Water for Cell Culture (WFI)

Size: 10 Liter 2D Bag

20 Liter 2D Bag

Product Description:

1A Clean Water for Injection (WFI) Quality Water and 1A Clean Water for Cell Culture (WFI) Quality Water are prepared by ultrafiltration, reverse osmosis, electro-de ionization and sterile filtration from a Water for Injection (WFI) system that meets or exceeds current EP/USP specifications.

Storage conditions:

Store at room temperature

Suitability / Areas of use:

Both products are suitable for compounding cell culture media and buffers, as well as a diluent. For research and manufacturing use only. These products are not intended for direct therapeutic use in humans or animals. Not for consumption.

Release Criteria:**Sterility:**

Pass – according to current EP – full product specific verification not included

Conductivity

≤ 1.1 µS/cm measured online in loop at 20 °C

Endotoxin

< 0.25 EU/ml

TOC

< 0.5 mg/L measured in loop

Nitrates

≤ 0.2 ppm